
DESCRIPTION
The Akreos™ Single Use Insertion Device is used for folding and delivering the Bausch & Lomb one-piece Akreos™ Advanced Optics Aspheric Lens into the eye. The device provides a small tubular pathway in which the Akreos Advanced Optics Aspheric Lens can be placed into the eye with one continuous forward motion. The Akreos Single Use Insertion Device consists of two parts:

1. Syringe shaped tube (includes lens loading deck and lens cap) and plunger (Fig. 1).

2. Transition cell (Fig. 2).

The Akreos Single Use Insertion Device is a sterile (EO), disposable plastic device, designed for single use only.

INDICATIONS

The Akreos Single Use Insertion Device should only be used with the Bausch & Lomb Akreos Advanced Optics Aspheric Lens from 0 to 30 diopters in power. Contact Bausch & Lomb for additional qualified product.

INSTRUCTIONS FOR USE

1. Open the package using standard sterile procedures. Place contents onto the sterile field.
2. Fill the rectangular end of the transition cell with approximately 0.2 ml of a Bausch & Lomb viscoelastic, using the line as a guide (Fig. 3).

3. Select the syringe-shaped tube and plunger. Pull the plunger until the tip of the plunger is just visible in the lens loading deck. Place a drop of viscoelastic directly in the center of the loading deck. Open the vial and remove the lens holder. Position the lens holder with its jaws facing down (Fig. 4a). With a blunt, non-toothed forceps, remove the lens from the lens holder by pushing the lens down. Transfer the lens to the recessed area in the injector loading deck (Fig. 4a), ensuring that the lens is not flipped during the transfer. The lens is loaded in correct orientation if its haptic features are as shown in Fig. 4b.

4. Advance the plunger until the fork fully captures the lens edge. Place a drop of viscoelastic in the center of the lens. Verify that the lens is inside of the recessed area of the loading deck, prior to closing the lens cap.

5. Close the lens cap (Fig. 5).

6. Align the key on the lens cap with the slot in the transition cell (Fig. 6). Push together with a firm motion until the tube with the cap and transition cell snap together (Fig. 7). Viscoelastic material will flow and fill both the lens chamber and the transition cell. Once assembled, pull on the transition cell to verify that it cannot be separated from the tube.

7. Fill the distal end of the device with viscoelastic material or balanced salt solution to reduce the possibility of introducing air pockets into the eye during implantation.

8. The lens is now ready to be delivered. Turn the device over, so that the bevel is face down, and insert the beveled portion of the transition cell through the incision in the eye. When the tip is positioned, apply continuous uniform forward motion to the plunger. A slow smooth forward motion of the plunger is essential for proper transition of the lens into position. Reverse movement of the plunger may result in the loss of the control of the lens or damage to the lens. Continue forward motion until the lens is fully expressed from the tip.

9. The edge of the unfolded lens may remain engaged in the plunger fork until it is released into the eye. To disengage the lens, release pressure on the plunger end allowing it to retract into the tip.

10. Discard the device after use.

PRECAUTIONS

Store the Akreos Single Use Insertion Device sterile sealed package at or below 30°C/86°F. Contents are sterile unless package is opened or damaged. Do not resterilize. Do not reuse. Single use only. Discontinue lens delivery if excessive force is needed to advance the plunger and lens.

MEDICAL DEVICE REUSE STATEMENT

If this product is reprocessed and/or re-used, Bausch + Lomb cannot guarantee the functionality, material structure, or cleanliness or sterility of the product. Re-use could lead to illness, infection and/or injury to the patient, user and, in extreme incidents, death. This product is labeled as 'single use' which is defined as a device intended to be used once only for a single patient.

CAUTION: Federal Law (USA) restricts this device to sale by, or on the order of, a physician.**SURGICAL PROCEDURE**

Preoperative surgical technique is the responsibility of the individual surgeon. Surgeons must determine the suitability of any particular procedure based upon their medical training and experience.

WARRANTY AND LIMITATIONS OF LIABILITY

Bausch & Lomb Incorporated warrants that the Akreos Single Use Insertion Device, when delivered, will conform to the manufacturer's then current version of the published specifications for such Akreos Single Use Insertion Device in all material respects and shall be free from defects in material or workmanship for a period equal to the Akreos Single Use Insertion Device expiration date. Bausch & Lomb Incorporated excludes all other warranties, whether expressed, implied or by operation of law or otherwise, including but not limited to any implied warranties of merchantability or fitness for a particular purpose. Bausch & Lomb Incorporated shall not be liable for any incidental, consequential or exemplary loss, damage or expense, directly or indirectly arising from the use of the Akreos Single Use Insertion Device even if Bausch & Lomb Incorporated has been advised of the possibility of such loss, damage or expense.

RETURN GOODS POLICY

All products returned to Bausch & Lomb Incorporated must be accompanied by a Return Goods Authorization Number. Call your customer service affiliate to obtain a Return Goods Authorization Number.

ADVERSE REACTION REPORTING

Adverse reactions and/or potentially sight-threatening complications that may be reasonably regarded as Akreos Single Use Insertion Device related and that were not previously expected in nature, severity or degree of incidence should be reported to your customer service affiliate.

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