Bausch & Lomb SofPort

Easy-Load Lens Delivery System

EN

DESCRIPTION

The EASY-LOAD LENS DELIVERY SYSTEM is a device used for folding and delivering the LI61AO, LI61AOV, or LI61SE 3-piece lens into the eye. The EASY-LOAD LENS DELIVERY SYSTEM consists of a syringe shaped body (1-A) and tip with a plunger (1-B), drawer (1-C), and haptic puller (1-D) shown in FIG. 1. The EASY-LOAD LENS DELIVERY SYSTEM is a sterile, disposable plastic device, designed for single use only. Combined with the lens case, FIG. 2, it provides fixture loading with no lens contact for the LI61AO and LI61AOV lens models or the manual load of the LI61SE lens model per FIG. 4, with a small tubular pathway in which the lens can be placed into the eye with one continuous forward motion.

INDICATIONS

The Easy-Load Lens Delivery System is indicated for the folding and injection of Bausch + Lomb intraocular lenses identifying the Easy-Load Lens Delivery System in their approved labeling.

INSTRUCTIONS FOR USE

1. Open the package using standard sterile procedures. Place the contents onto the sterile field.
2. Confirm that the haptic puller is fully installed into the body.
3. Confirm that the plunger is in the start position shown in FIG. 1.
4. While keeping the device level, apply Bausch & Lomb viscoelastic (3-A) through the hole of the middle drawer finger and to the floor of the loading area extending under the lens track edge as shown in FIG. 3.
5. For the LI61SE, manually place the lens in the EZ-28 loading area using non-serrated forceps. Position the end of the trailing haptic over the plunger tip and the leading haptic under the lens track edge and over the haptic puller, as shown in FIG. 4. Now proceed directly to Step 10 to complete the LI61SE lens loading.
6. For the LI61AO & AOV, separate the lens retainer (5-A) from the lens cover (5-B) with a pivot and pull motion as shown in FIG. 5.
7. Place the lens retainer in the load area. Engage the lens retainer to the body with a slow pivoting motion as shown in FIG. 6.
8. Actuate the drawer by pressing up (7-A) on the drawer stop arm as shown in FIG. 7. Push slowly forward with the thumb until the lens retainer has stopped (8-A) the drawer movement as shown in FIG. 8. The lens is now ready to be removed from the lens retainer.
9. Pivot the lens retainer slightly away from drawer stop and lift away from the body as shown in FIG. 9. Discard the lens retainer after removal.
10. Slowly close the drawer until the snap closure mechanism has engaged as shown in FIG. 10. **NOTE:** Closing the drawer compresses the lens for delivery. Do not close the drawer, thereby compressing the lens, until immediately before insertion. Closure of the drawer can be verified by the inability to pull the drawer back from its snapped position.
11. Advance the plunger forward, the haptic puller will begin to move away from the body tip. **NOTE:** WHILE ADVANCING THE PLUNGER DO NOT RETRACT UNTIL THE OPTIC IS DELIVERED AND FREE OF THE TIP. As the plunger moves forward it will stop at a specific position (11-A) as shown in FIG. 11. Straighten the leading haptic by removing the haptic puller. The straightened haptic will be inside the body tip. Discard the haptic puller.
12. Fill the distal end of the EASY-LOAD LENS DELIVERY SYSTEM with viscoelastic material or balanced salt solution to reduce the possibility of introducing air pockets into the eye during lens implantation.
13. Insert the beveled portion of the tip through the incision in the eye. When the tip is positioned, apply uniform forward pressure on the plunger until the lens is fully expressed from the tip and delivered into the capsular bag. Once pressure is released from the plunger, it will retract automatically to engage the trailing haptic. When the trailing haptic is engaged, uniform forward pressure is once again applied on the plunger until the trailing haptic is completely delivered into the capsular bag.
14. Discard the EASY-LOAD LENS DELIVERY SYSTEM after use.

PRECAUTIONS AND WARNINGS

Store the EASY-LOAD LENS DELIVERY SYSTEM's sterile sealed package at room temperature (<30°C/86°F). Contents are sterile unless package is opened or damaged. Do not resterilize. Folding and compression of the lens should be done just prior to insertion and delivery. Failure to follow the Instructions for Use may result in patient injury.

SURGICAL PROCEDURE

Proper surgical technique is the responsibility of the individual surgeon. Surgeons must determine the suitability of any particular procedure based upon their medical training and experience.

WARRANTY AND LIMITATIONS OF LIABILITY

Bausch & Lomb Incorporated warrants that the EASY-LOAD LENS DELIVERY SYSTEM, when delivered, will conform to the manufacturer's then current version of the published specifications for such EASY-LOAD LENS DELIVERY SYSTEM in all material respects and shall be free from defects in material or workmanship for a period equal to the EASY-LOAD LENS DELIVERY SYSTEM'S expiration date.

BAUSCH & LOMB INCORPORATED EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED, IMPLIED OR BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS. BAUSCH & LOMB INCORPORATED SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL OR EXEMPLARY LOSS, DAMAGE OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM THE USE OF THE EASY-LOAD LENS DELIVERY SYSTEM EVEN IF BAUSCH & LOMB INCORPORATED HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, DAMAGE OR EXPENSE.

RETURN GOODS POLICY

All products returned to Bausch & Lomb Incorporated must be accompanied by a Return Authorization number. Call 1-800-338-2020 for a return authorization and full policy information.

FIG. 1


FIG. 2


FIG. 3


FIG. 4


FIG. 5


FIG. 6


FIG. 7


FIG. 8


FIG. 9


FIG. 10



FIG. 11

Bausch & Lomb Incorporated
Rochester, NY 14609 USA

EC REP Bausch & Lomb Incorporated
106 London Road
Kingston-upon-Thames, KT2 6TN, UK

Manufacturing Site:
Bausch & Lomb Incorporated
21 Park Place Blvd. North
Clearwater, FL 33759 USA

®/TM are trademarks of Bausch & Lomb Incorporated or its affiliates.
© Bausch & Lomb Incorporated. All rights reserved.
U.S. Patents: 5944725, 6491697 and 7422604. Other patents pending.

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*TM are trademarks of Bausch + Lomb Incorporated or its affiliates.
© Bausch + Lomb Incorporated. All rights reserved.