CRYSTALCERT™ DELIVERY SYSTEM

DESCRIPTION

The Crystalsert™ delivery system is a device used for folding and delivering the Crystallens accommodating intraocular lens (AT-52SE, AT-50SE, HD520, HD500, AT50AO, AT52AO), Trulign™ Toric, and other IOLs indicating use of the Crystalsert delivery system in their approved labeling into the eye. The Crystalsert delivery system consists of a syringe shaped body (1A) and tip (1B) with a plunger (1C) and drawer (1D) shown in FIG. 1. The Crystalsert delivery system is a sterile, disposable plastic device, with a small tubular pathway in which the lens can be placed into the eye with one continuous motion, designed for single use only.

INDICATIONS

The Crystalsert™ delivery system is intended to be used to fold and deliver the Crystallens accommodating intraocular lens and other intraocular lens identifying the Crystalsert delivery system in their approved labeling.

INSTRUCTIONS FOR USE

1. Open the package using standard sterile procedures. Place the contents onto the sterile field.
2. Confirm that the plunger is in the start position (1E) as shown in FIG. 1.
3. While keeping the device level, apply Bausch & Lomb Armvisc Plus viscoelastic to the floor of the loading area (2A) and the area extending under the lens track edge (2B) as shown in FIG. 2.
4. Place the lens in the load area as shown in FIG. 3, with optic in contact with the body floor (3A). Note that the leading circular haptic (3B) is under the lens track edge and that the leading plate haptic is tangent to the body edge (3C). The trailing haptic should be in contact with the lens track edge (3D). Ensure that the plate haptics (3E) are also in contact with the body floor (3A). The trailing circular haptic lobe (3F) must be positioned as shown for correct lens orientation.
5. Actuate the drawer by pressing up on the drawer stop arm (4A) as shown in FIG. 4.
6. Slowly close the drawer until the snap closure mechanism has engaged as shown in FIG. 5.
7. **NOTE:** Closing the drawer compresses the lens for delivery. Do not close the drawer, thereby compressing the lens, until immediately before insertion. A properly closed drawer cannot be pulled back from its snapped position. The trailing haptic (6A) should clearly be seen in the opening shown in FIG. 6.
8. Advancing the plunger forward moves the lens towards the body tip. As the plunger moves approximately .75 inches or 19 mm forward it will stop at a detent position (7A) as shown in FIG. 7.
9. At the detent position fill the distal end of the Crystalsert™ delivery system with Bausch & Lomb Armvisc Plus viscoelastic or balanced salt solution to reduce the possibility of introducing air pockets into the eye during lens implantation.
10. Insert the beveled portion of the tip through the incision in the eye, with the bevel facing down. When the tip is positioned, apply uniform forward pressure on the plunger until the lens is fully expressed from the tip and delivered into the capsular bag.
11. Discard the Crystalsert delivery system after use.

PRECAUTIONS AND WARNINGS

Store the Crystalsert™ delivery system's sterile sealed package at room temperature (<30 °C/86 °F). Contents are sterile unless package is opened or damaged. Do not resterilize. Folding and compression of the lens should be done just prior to insertion and delivery. Failure to follow the Instructions for Use may result in patient injury.

SURGICAL PROCEDURE

Proper surgical technique is the responsibility of the individual surgeon. Surgeons must determine the suitability of any particular procedure based upon their medical training and experience.

WARRANTY AND LIMITATIONS OF LIABILITY








Bausch & Lomb Incorporated warrants that the Crystalsert™ delivery system, when delivered, will conform to the manufacturer's then current version of the published specifications for such Crystalsert™ delivery system in all material respects and shall be free from defects in material or workmanship for a period equal to the Crystalsert™ delivery system's expiration date.


BAUSCH & LOMB INCORPORATED EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED, IMPLIED OR BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS. BAUSCH & LOMB INCORPORATED SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL OR EXEMPLARY LOSS, DAMAGE OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM THE USE OF THE CRYSTALCERT™ DELIVERY SYSTEM EVEN IF BAUSCH & LOMB INCORPORATED HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSS, DAMAGE OR EXPENSE.

RETURN GOODS POLICY

All products returned to Bausch & Lomb Incorporated must be accompanied by a Return Authorization number. Call 1-800-338-2020 for a return authorization and full policy information.

*/™ are trademarks of Bausch & Lomb Incorporated or its affiliates.
© Bausch & Lomb Incorporated. All rights reserved.
U.S. Patent: 5944725. Other patents pending.

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