**Crystalsert®**

Delivery System

DISPOSITIF DE PLACEMENT CRYSTALSET*

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6. Assurez-vous que la plante forward moves the lens towards the body tip. As the plunger moves approximately .75 inches or 19 mm forward it will stop at a detent position (7A) as shown in Figure 6.

7. NOTE: Closing the closure compresses the lens for delivery. Do not close the closure, thereby compressing the lens, until immediately before insertion. A properly closed closure cannot be pulled back from its snapped position. The trailing haptic (6A) should close when the closure is open.

8. Après avoir fermé la plante forward moves the lens towards the body tip. As the plunger moves approximately .75 inches or 19 mm forward it will stop at a detent position (7A) as shown in Figure 6.

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10. Insert the bevelled portion of the tip through the incision in the eye, with the bevel facing down. When the tip is positioned, apply uniform forward pressure on the plunger until the lens is fully expressed from the tip and delivered into the capsular bag.

11. Discard the Crystalsert delivery system after use.

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MEDICAL DEVICE RE-USE STATEMENT
If this product is reproduced and/or reused, Bausch + Lomb cannot guarantee that its functionality, material, sterility or cleanliness of the product. Re-use could lead to illness, infection and/or injury to the patient or user, and, in extreme incidents, death. It is the responsibility of any particular procedure based upon their medical training and expertise.

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WARRANTY AND LIMITATIONS OF LIABILITY
Bausch + Lomb Incorporated warrants that the Crystalsert® delivery system, when delivered, will conform to the manufacturer's then current version of the published specifications for such Crystalsert® delivery system in all material respects and shall be free from defects in materials or workmanship for a period equal to the Crystalsert® delivery system's expiration date.

BAUSCH & LOMB INCORPORATED EXCLUSIONS OF WARRANTIES
EXPRESSION, IMPLIED OR BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO THE EXCLUSIONS OF LIABILITY FOR FAULT, DAMAGE OR INCONVENIENT SHALL NOT BE APPLIED FOR ANY INCIDENTAL, CONSEQUENTIAL OR EXEMPLARY LOSS, DAMAGE OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM THE USE OF THE CRYSTALSET® DELIVERY SYSTEM EVEN IF BAUSCH & LOMB INCORPORATED HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH LOSSES.

RETURNS POLICY
All products returned to Bausch + Lomb Incorporated must be accompanied by a Return Authorization number. Call 1-800-338-2020 for a return authorization and full policy information.

PROCEDURE FOR RETURNS
Bausch + Lomb Incorporated warrants that the Crystalsert® delivery system, when delivered, will conform to the manufacturer's then current version of the published specifications for such Crystalsert® delivery system in all material respects and shall be free from defects in materials or workmanship for a period equal to the Crystalsert® delivery system's expiration date.

CRYSTALSET® DELIVERY SYSTEM
Bausch + Lomb Incorporated warrants that the Crystalsert® delivery system, when delivered, will conform to the manufacturer's then current version of the published specifications for such Crystalsert® delivery system in all material respects and shall be free from defects in materials or workmanship for a period equal to the Crystalsert® delivery system's expiration date.

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